



Alcoholic rub for hygienic and surgical hand disinfection.

desderman® pure gel

Our Plus

- meets the FDA clinical efficacy standards (TFM)
- tested in accordance with European standards
- suitable for hygienic and surgical hand disinfection
- colour and fragrance free, therefore hypoallergenic
- excellent skin compatibility and pleasant skin feeling due to an established system of refatting agents

Application areas

- hygienic hand disinfection (EN 1500): Rub 3 ml into clean dry hands, contact time: 30 sec.
- surgical hand disinfection (EN 12791): Rub a sufficient amount into clean dry hands and forearms, contact time 90 sec.

Instructions for use

desderman® pure gel is to be used undiluted. During the contact time for hygienic and surgical hand disinfection, desderman® pure gel should be applied in portions so that the hands respectively the hands and forearms remain well moistened during the entire application time.

Microbiological efficacy

Efficacy	Concentration	Contact time
bactericidal	ready-to-use	30 sec.
Escherichia coli	ready-to-use	30 sec.
Klebsiella pneumoniae	ready-to-use	30 sec.
MRSA EN13727	ready-to-use	30 sec.
Staphylococcus aureus	ready-to-use	30 sec.
yeastcidal EN13624	ready-to-use	15 sec.

Efficacy	Concentration	Contact time
fungicidal EN13624	ready-to-use	60 sec.
virucidal against enveloped viruses in accordance with DVV (German Registered Association for Combating Viral Diseases)//RKI Guideline	ready-to-use	30 sec.
Vaccinia virus	ready-to-use	30 sec.
limited spectrum virucidal activity EN14476	ready-to-use	2 min.
virucidal EN14476	ready-to-use	2 min.
Adenovirus EN14476	ready-to-use	2 min.
Norovirus EN14476	ready-to-use	15 sec.
Poliovirus EN14476	ready-to-use	90 sec.
Rotavirus EN14476	ready-to-use	15 sec.

Application area	Concentration	Contact time
hygienic hand disinfection EN1500	ready-to-use	30 sec.
surgical hand disinfection EN12791	ready-to-use	90 sec.

Certificates

- VAH certificate



desderman[®] pure gel

Product data

Composition

100 g solution contains the following active substances: Active ingredients: 78.2 g ethanol 96 %, 0.1 g biphenyl-2-ol.

Chemical-physical data

Color	colourless
Density	approx. 0,83 g/cm ³ / 20 °C
Flash point	13 °C / Method : DIN 53213, Part 1
Form	viscous
pH	Not applicable
Viscosity, dynamic	700 - 1.300 mPa*s / 20 °C / Method : DIN 53019

Special advice

Use disinfectants safely. Always read the label and product information before use.

Do not use after expiry date.

Information for order

Item	Delivery form	Item no.
desderman pure gel TLD 1 l FL	6/Carton	on request
desderman pure gel INT 100 ml FL	30/Carton	on request
desderman pure gel 500 ml FL	20/Carton	on request
desderman pure gel 1 l FL	10/Carton	on request
desderman pure gel -INT- 1 l EF	10/Carton	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

Application aids

Application aids	Item no.
Dispenser KHK 1000 (approx. 0.75-1.5ml per stroke)	669700
Dispenser KHK 500 (approx. 0.75-1.5ml per stroke)	669600
Product dispenser long arm 1000 ml	669710
Product dispenser long arm 500 ml	669610
schülke touchless D 1	668300
sm 2 500 (approx- 1.0-3.0ml per stroke)	668600
sm 2 Universal (approx. 1.0 - 3.0 ml per stroke)	668500

Related Products

- desderman[®] pure
- sensiva[®] dry skin balm
- sensiva[®] protective cream
- sensiva[®] protective emulsion
- sensiva[®] regeneration cream
- sensiva[®] wash lotion

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: www.schuelke.com.

For individual questions: Customer Sales Service Phone: +49 40 52100-666E-Mail: info@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com